

Carboxylierte Styrol-Butadien-Latizes	Gehalt Styrol (%)	Brookfield-Viskosität (20 rpm, 25°C, mPa.s)
EUOPRENE® Latice 405	40	300
EUOPRENE® Latice 406	40	300
EUOPRENE® Latice 440	60	300
EUOPRENE® Latice 455	47	300
EUOPRENE® Latice 5583	40	300
EUOPRENE® Latice 5584	60	400
EUOPRENE® Latice 5585	47	300
EUOPRENE® Latice 5587	75	600
EUOPRENE® Latice 5588	50	350
EUOPRENE® Latice 8435	69	600
EUOPRENE® Latice 8487	67	500